



Head Office

1406 Coventry Road, South Yardley
Birmingham, B25 8 AE .UK
Tel: +44 121 765 0214
Fax: +44 121 765 0215
COMP. REG. No.: 4917976

Regional Office

P.O.Box: 9442
Saif Zone, Sharjah - U.A.E.
Tel.: 06-5578611
Fax.: 06-5578612

www.halleymedical.co.uk

Email: sales@halleymedical-uk.ae

MEDICAL PRODUCTS 2011

MEDICAL PRODUCTS



Halley Medical Supplies Ltd.

t h e g e n t l e t o u c h

2011

HALLEY MEDICAL SUPPLIES LTD. UK

- Medical Consumables
- Wound Care Products
- Disposable Gloves
- Medical disposables
- Syringes & Needles
- Surgical Gowns & Drapes
- Surgical Face Masks & Caps



CONTENTS

MEDICAL CONSUMABLES

Adhesive Tapes	04-07
Bandages	08-14
Disposable Gloves	15-21
Cotton Products	22-25
Dressings	26-27
Gauze Products	28-30
Assorted Products	31-33
Diagnostic Products	34-38
First Aid Boxes	39-41

MEDICAL DISPOSABLES

Syringes & Needles	44-46
I.V Catheters & Sets	47-52
Disposal Containers	53
Special Kits	54-56
Assorted Products	57-65
Masks & Caps	66-71
Gauze products	72-74
Gowns & Drapes	75-76



PREFACE

*H*ealthcare requirements and needs in the entire World have been always increasingly at advanced rate; hence Halley Medical is committed to provide the best healthcare solutions. Our mission is to service the Human being well and with all relevant needs, we certainly assure you of highest and unique quality, as continuous product development and enhancement has given HMS the edge over the competition, beside the production flexibility to meet all your needs, as well as our initiative and creative ideas are endless. We focus strategically on the remarkable growth Worldwide by continuously increasing the HMS vast range of products, which includes: Medical disposables & consumables like adhesive tapes, bandages, wound plasters, cotton products, gauze products, gloves, diagnostic products, face masks, head caps, surgical gowns & drapes, and lap sponges, beside surgical instruments, laboratory consumables, hospital textiles and uniforms, and rehabilitation products.

ZINC OXIDE PLASTER

ZINC OXIDE ADHESIVE TAPE

Surgical adhesive tape, made of pure cotton fabric,
 Skin coloured with zig-zag edges for easy tear off,
 With reliable adhesive property and high breathability.



Consumer Pack

Consumer Pack

Cat.No	Size	Pcs/Pkt	Pcs/Carton
28.0100	1.25 cm x 5 y	24	480
28.0101	2.5 cm x 5 y	12	360
28.0102	5 cm x 5 y	6	180
28.0103	7.5 cm x 5 y	6	180
28.0104	10 cm x 5 y	6	144
28.0110	1.25 cm x 10 y	24	480
28.0111	2.5 cm x 10 y	12	360
28.0112	5 cm x 10 y	6	180
28.0113	7.5 cm x 10 y	6	180
28.0114	10 cm x 10 y	6	144

NON-WOVEN TAPE

NON-WOVEN SURGICAL TAPE

Surgical adhesive tape, made of hypoallergenic
 non-woven material for sensitive skin,
 High permeability to the air, easy to tear off,
 Skin friendly; painless removal.



Medical Pack

Consumer Pack

Cat.No	Size	Pcs/Pkt	Pcs/Carton
28.0200	1.25 cm x 5 y	24	480
28.0201	2.5 cm x 5 y	12	360
28.0202	5 cm x 5 y	6	180
28.0210	1.25 cm x 10 y	24	480
28.0211	2.5 cm x 10 y	12	360
28.0212	5 cm x 10 y	6	180

SILK ADHESIVE TAPE

SILK ADHESIVE TAPE

Hypoallergenic surgical adhesive tape for sensitive skin,
Highly air permeable; easy to tear off,
Can be removed painlessly and without leaving residues.



Medical Pack


Cat.No	Size	Pcs/Pkt	Pcs/Carton
28.0300	1.25 cm x 5 y	24	480
28.0301	2.5 cm x 5 y	12	360
28.0302	5 cm x 5 y	6	180
28.0310	1.25 cm x 10 y	24	480
28.0311	2.5 cm x 10 y	12	360
28.0312	5 cm x 10 y	6	180

Consumer Pack

P.E. TRANSPARENT SURGICAL TAPE

P.E. ADHESIVE TAPE

Hypoallergenic surgical transparent adhesive tape,
High permeability to the air and water vapour,
Easy to tear off, with reliable adhesive property.



Medical Pack

Cat.No	Size	Pcs/Pkt	Pcs/Carton
28.0400	1.25 cm x 5 y	24	480
28.0401	2.5 cm x 5 y	12	360
28.0402	5 cm x 5 y	6	180
28.0410	1.25 cm x 10 y	24	480
28.0411	2.5 cm x 10 y	12	360
28.0412	5 cm x 10 y	6	180

Consumer Pack

RUBBER ELASTIC BANDAGE

Multipurpose elastic bandage,
 reinforced with natural rubber for
 high stretch compression property,
 Skin colour, washable and reusable,
 Individually wrapped, with 2 metal clips,



Cat.No	Size	Pcs/Pkt	Pcs/Carton
29.0102	2" x 4.5 m	6	360
29.0103	3" x 4.5 m	6	240
29.0104	4" x 4.5 m	6	180
29.0105	5" x 4.5 m	6	162
29.0106	6" x 4.5 m	6	132

COTTON CREPE ELASTIC BANDAGE

Excellent stretch, allow free movement,
 Comfortable, conforms well to all body area,
 White colour, washable and reusable,
 Individually wrapped, with 2 metal clips.



Cat.No	Size	Pcs/Pkt	Pcs/Carton
29.0202	5 cm x 4.5 m	6	360
29.0203	7.5 cm x 4.5 m	6	240
29.0204	10 cm x 4.5 m	6	180
29.0205	12.5 cm x 4.5 m	6	162
29.0206	15 cm x 4.5 m	6	132

COTTON ADHESIVE ELASTIC BANDAGE

Adhesive compression bandage,
 Latex Free, No allergic reaction,
 Soft, flexible, breathable and comfortable,
 Available colours: White & Skin coloured
 individually wrapped.



	Cat.No	Size	Pcs/Pkt	Pcs/Carton
Skin coloured	29.0303	2" x 4.5 m	6	216
	29.0304	3" x 4.5 m	6	162
	29.0305	4" x 4.5 m	6	108
White colour	29.0356	2" x 4.5 m	6	216
	29.0357	3" x 4.5 m	6	162
	29.0358	4" x 4.5 m	6	108

ELASTIC COHESIVE BANDAGE

Consist of poly cotton fabrics,
 Cohesive properties (Self-adhesive) allowing
 painless removal,
 Breathable and Comfortable,
 Excellent stretch, allow free movement,
 individually wrapped.



	Cat.No	Size	Pcs/Pkt	Pcs/Carton
Skin coloured	29.0403	5 cm x 4.5 m	6	216
	29.0404	7.5 cm x 4.5 m	6	162
	29.0405	10 cm x 4.5 m	6	108
White colour	29.0464	5 cm x 4.5 m	6	216
	29.0465	7.5 cm x 4.5 m	6	162
	29.0466	10 cm x 4.5 m	6	108

TUBULAR ELASTIC BANDAGE


Elastic and comfortable
easy to cut, and can be used on all parts of the body
effective supports for sprains, strains,,,,etc
provided in dispensing box.



Cat.No	Size	Pcs/Pkt	Pcs/Carton
29.0831	4.5 cm x 10 m	1 Roll	28
29.0833	6.25 cm x 10 m	1 Roll	20
29.0835	6.75 cm x 10 m	1 Roll	20
29.0837	7.5 cm x 10 m	1 Roll	20
29.0839	8.75 cm x 10 m	1 Roll	16
29.0841	10 cm x 10 m	1 Roll	16
29.0843	12 cm x 10 m	1 Roll	16

ELASTIC TUBULAR NET BANDAGE

High elasticity and breathability
dressing secure with free movement
Comfortable, conforms well to all body area,
provided in dispensing box.



Cat.No	Size	Pcs/Pkt	Pcs/Carton
29.0871	0.7 cm x 25 m (Size 0.5#)	1 Roll	80
29.0872	1.7 cm x 25 m (Size 1.0#)	1 Roll	60
29.0873	2.0 cm x 25 m (Size 2.0#)	1 Roll	60
29.0874	2.3 cm x 25 m (Size 3.0#)	1 Roll	42
29.0875	2.7 cm x 25 m (Size 4.0#)	1 Roll	42
29.0876	3.0 cm x 25 m (Size 5.0#)	1 Roll	36
29.0877	3.5 cm x 25 m (Size 5.5#)	1 Roll	36
29.0878	4.2 cm x 25 m (Size 5.8#)	1 Roll	36
29.0879	5.8 cm x 25 m (Size 6.0#)	1 Roll	32
29.0880	7.2 cm x 25 m (Size 7.0#)	1 Roll	32
29.0881	8.5 cm x 25 m (Size 8.0#)	1 Roll	32

CONFIRMING BANDAGE

Breathable and Comfortable,
 Bleached colour, Poly cotton material
 High elastic fixation bandage for dressing secure
 Individually wrapped,



Cat.No	Size	Pcs/Pkt	Pcs/Carton
29.0672	5 cm x 4.5 m	12	480
29.0674	7.5 cm x 4.5 m	12	324
29.0676	10 cm x 4.5 m	12	240
29.0678	15 cm x 4.5 m	12	180

SURGICAL STERILE GLOVES

Powdered & Powder Free

Natural Rubber Latex
 Micro Rough Surface for optimal instruments grip,
 Contain low protein allergens for sensitive skin,
 Dusted with absorbable Corn-Starch Powder (U.S.P)
 Manufacturer in accordance with BS 4005/ EN 455/ ASTM D-3577
 Packed in pairs in sterile peel pouch
 Beaded Cuff, Single use only.



Cat.No (Powdered)	Cat.No (Powder Free)	Size	Pcs/Pkt	Pcs/Carton
21.9551	21.9571	6.0	100 (50 Pairs)	800
21.9552	21.9572	6.5	100 (50 Pairs)	800
21.9553	21.9573	7.0	100 (50 Pairs)	800
21.9554	21.9574	7.5	100 (50 Pairs)	800
21.9555	21.9575	8.0	100 (50 Pairs)	800
21.9556	21.9576	8.5	100 (50 Pairs)	800

LATEX EXAMINATION GLOVES

Lightly Powdered

Natural Rubber Latex
 Dusted with absorbable Corn-Starch Powder (U.S.P)
 Contain low protein allergens,
 Meet ASTM D 3578 and fulfils to CE standards
 Beaded Cuff,
 Smooth,
 Ambidextrous
 Non-Sterile,
 Single use only.




Cat.No	Size	Pcs/Pkt	Pcs/Carton
21.0101	X-Small	100	1000
21.0102	Small	100	1000
21.0103	Medium	100	1000
21.0104	Large	100	1000
21.0105	X-Large	100	1000

LATEX EXAMINATION GLOVES

Powder Free

For Sensitive skin
 Natural Rubber Latex
 Contain low protein allergens,
 Meet ASTM D 3578 and fulfils to CE standards
 Beaded Cuff,
 Textured surface for improved grip,
 Ambidextrous
 Non-Sterile,
 Single use only.




Cat.No	Size	Pcs/Pkt	Pcs/Carton
21.0121	X-Small	100	1000
21.0122	Small	100	1000
21.0123	Medium	100	1000
21.0124	Large	100	1000
21.0125	X-Large	100	1000



DISPOSABLE VINYL GLOVES

Lightly Powdered

Gloves are lightly powdered with absorbable corn starch (U.S.P),
 Manufacturer in accordance with EN 455/ ASTM D-5250
 Gloves are made of 100% DINP material
 DINP is a high quality, environmentally friendly, and is a low toxicity plasticizers,
 Gloves are 100% DEHP Free,
 Colors: Clear – Blue
 Smooth, Ambidextrous, Non-Sterile, Single use only.




Cat.No (Clear color)	Size	Pcs/Pkt	Pcs/Carton
21.0211	X-Small	100	1000
21.0212	Small	100	1000
21.0213	Medium	100	1000
21.0214	Large	100	1000
21.0215	X-Large	100	1000

Cat.No (Blue color)	Size	Pcs/Pkt	Pcs/Carton
21.0231	X-Small	100	1000
21.0232	Small	100	1000
21.0233	Medium	100	1000
21.0234	Large	100	1000
21.0235	X-Large	100	1000

DISPOSABLE VINYL GLOVES

Powder Free

for those allergic to powder
 Manufacturer in accordance with EN 455/ ASTM D-5250
 Gloves are made of 100% DINP material
 DINP is a high quality, environmentally friendly, and is a low toxicity plasticizers,
 Gloves are 100% DEHP Free,
 Colors: Clear – Blue
 Ambidextrous, Non-Sterile, Single use only.



Cat.No (Clear color)	Size	Pcs/Pkt	Pcs/Carton
21.0216	X-Small	100	1000
21.0217	Small	100	1000
21.0218	Medium	100	1000
21.0219	Large	100	1000
21.0220	X-Large	100	1000

Cat.No (Blue color)	Size	Pcs/Pkt	Pcs/Carton
21.0236	X-Small	100	1000
21.0237	Small	100	1000
21.0238	Medium	100	1000
21.0239	Large	100	1000
21.0240	X-Large	100	1000

POLYTHENE GLOVES

PE Gloves

Embossed for better grip

available in PE bag and in carton box

Ambidextrous

Non-Sterile,

Single use only.



Cat.No LDPE/bag	Cat.No LDPE/Box	Size	Pcs/Pkt	Pcs/Carton
21.0332	21.0337	Small	100	10,000
21.0333	21.0338	Medium	100	10,000
21.0334	21.0339	Large	100	10,000
Cat.No HDPE/bag	Cat.No HDPE/Box	Size	Pcs/Pkt	Pcs/Carton
21.0342	21.0347	Small	100	10,000
21.0343	21.0348	Medium	100	10,000
21.0344	21.0349	Large	100	10,000

NITRILE EXAMINATION GLOVES

LATEX FREE & POWDER FREE

Ideal for those with Latex Allergies

Manufacturer in accordance with

EN 455/ ASTM D-6319

available colours: white and blue

Textured surface for improved grip,

Ambidextrous, Non-Sterile, Single use only.



White color-Powder Free	Cat.No	Size	Pcs/Pkt	Pcs/Carton
	21.0416	X-Small	100	1000
	21.0417	Small	100	1000
	21.0418	Medium	100	1000
	21.0419	Large	100	1000
	21.0420	X-Large	100	1000
Blue color-Powder Free	Cat.No	Size	Pcs/Pkt	Pcs/Carton
	21.0436	X-Small	100	1000
	21.0437	Small	100	1000
	21.0438	Medium	100	1000
	21.0439	Large	100	1000
	21.0440	X-Large	100	1000

COTTON ROLLS

Absorbent Cotton Wool

100 % Cotton Wool,

Hygienic purity with high absorption capacity



Cat.No	weight/Roll	Rolls/Carton
30.0105	25 grs	400
30.0107	50 grs	200
30.0109	100 grs	200
30.0115	300 grs	50
30.0117	500 grs	30
30.0119	1000 grs	15

COTTON WOOL ZIG ZAG

Easy to tear , absorbent Cotton Wool

100 % Cotton Wool,

Hygienic purity with high absorption capacity

Packed in plastic bag,



Cat.No	weight/Pkt	Pkt/Carton
30.0227	25 grs	360
30.0229	50 grs	240
30.0233	100 grs	120

COTTON BALLS

100 % Cotton Wool,
 available colours: white and Coloured
 Hygienic purity with high absorption capacity
 Packed in plastic zip bag,



	Cat.No	weight/Pkt	Pkt/Carton
White	30.0503	40 grs	72
	30.0505	100 grs	24
Coloured	30.0507	40 grs	72
	30.0509	100 grs	24

COTTON BUDS

100 % Cotton tips and polypropylene sticks,
 Sticks colours: white and coloured
 Suitable for cosmetic use
 Packed in PVC can,



Cat.No	Pcs/Pkt	Pkt/Carton
30.0301	100	240
30.0302	200	240



WOUND ADHESIVE PLASTER

Breathable-Flexible with highly absorbent non-stick pad for gentle removal,
 Available materials : PVC , Fabric (Cotton) and Stretched Fabric,
 Available types: Strip Plaster and round Plaster



	Cat.No	Size	Pcs/Pkt	Pcs/Carton
72x19mm	25.9103	PVC Plaster	100	10,000
	25.9105	Fabric Plaster	100	10,000
	25.9107	Stretched Fabric Plaster	100	10,000
72x25mm	25.9113	PVC Plaster	100	10,000
	25.9121	PVC Round Spots	20	6000
22 mm	25.9125	PVC Round Spots	100	10,000

SURGICAL DRESSING

Adhesive Surgical Dressing with highly absorbent pad,
 lifts off gently,
 Breathable-Flexible and comfortable,
 Sterilized and Packed into individual paper pouches,



Cat.No	adhesive Sheet	Absorbent Pad	Pcs/Pkt	Pcs/Carton
25.9304	6 x 8 cm	3 x 4 cm	50	1200
25.9308	10 x 10 cm	5 x 5 cm	50	1200
25.9312	10 x 15 cm	5 x 10 cm	50	1200
25.9316	10 x 20 cm	4 x 14 cm	50	1200
25.9320	10 x 25 cm	5 x 19 cm	50	1200



GAUZE SWAB-NON STERILE

100 % Cotton gauze swabs,
Soft & breathable
folded edges,
40's / 19 x 15 mesh,
100 Pcs /Pkt



	Cat.No	Size	Pcs/Pkt	Pcs/Carton
8 Plies	27.0402	5 x 5 cm	100	14,400
	27.0403	7.5 x 7.5 cm	100	6400
	27.0404	10 x 10 cm	100	3600
12 Plies	27.0602	5 x 5 cm	100	14,400
	27.0603	7.5 x 7.5 cm	100	6400
	27.0604	10 x 10 cm	100	3600

W.O.W GAUZE BANDAGE

Consist of Pure absorbent cotton gauze as per USP/BP,
Soft, comfortable, and high breathability,
12 Pcs, each individually wrapped in paper pack,
and all packed in a pouch.
40's/19 x 15, 100 Doz/Ctn



Cat.No	Size	Pcs/Pkt	Pcs/Carton
27.0901	2.5 cm x 5 m	12	1200
27.0902	5 cm x 5 m	12	1200
27.0903	7.5 cm x 5 m	12	1200
27.0904	10 cm x 5 m	12	1200
27.0906	15 cm x 5 m	12	1200

GAUZE ROLL

Made of 100 % absorbent cotton gauze,
High absorbency and breathability,
40's/19 x 15, 4 Plies, 36"x 100 yds,



Cat.No	Size	Pcs/Pkt	Pcs/Carton
27.0022	36" x 100 Yds	1	20

TRIANGULAR BANDAGE,

100 % Cotton
Forearm and wrist support
1 Pec/Poly Bag -12 Pcs/master poly bags - 240 Pcs/Ctn



Cat.No	Size	Pcs/Pkt	Pcs/Carton
24.0254	36" x 36"x 51"	12	240

DUST FACE MASK

Provides a comfortable fit, easy breathing, light weight, clear vision, Latex and Fiber free. multipurpose mask recommended for general healthcare practice, household , and light duty tasks,



Cat.No	Pcs/Pkt	Pcs/Carton
24.0305	50	1000

PRE-INJECTION SWABS

Alcohol Prep Pad – for external use only
Contents: Saturated with 70% Isopropyl Alcohol,
For single use only.



Cat.No	Pcs/Pkt	Pcs/Carton
24.0101	100	10,000
24.0102	200	20,000

STANDARD STETHOSCOPE

STANDARD STAINLESS STEEL STETHOSCOPE

HMS Stainless Steel Stethoscope, has been carefully designed for Excellent Performance,

Chest Piece made from high grade S.S for Optimal performance and rugged durability,

Extra deep open bell (1 38/”) for detection of low frequency sounds,

Ultra sensitive diaphragm (1 18/”) for greater amplification and crisper high frequency transmission,

Complete with tow pairs of PVC Ear Tips and with ID Patch.



Cat.No	Size	Pcs/Pkt	Pcs/Carton
22.0123	Adult	1	20
22.0152	Paediatric	1	20

MERCURIAL SPHYGMOMANOMETER

Desk Type Mercurial Sphygmomanometer

Cotton cuff, Grey colour case

D-ring cuff

0300- mmHg



Cat.No	Pcs/Pkt	Pkt/Carton
22.0237	1	20

ANEROID SPHYGMOMANOMETER

With Single Head Stethoscope
 Cotton grey colour cuff with D-Ring,
 300 mm Hg gauge



Cat.No	Pcs/Pkt	Pkt/Carton
22.0317	1	20

PEN LIGHTS

Durable ABS Plastic,
 With on/off button,
 Light weight,
 With pocket clip,
 Supplied with 2 x AAA batteries,



CLINICAL THERMOMETER

Oral , Rectal and Armpite types
 Prismatic shape (Oral & Rectal)
 Flat Shape (Armpite)
 C&F,
 Packed in plastic case



Cat.No	Type	Pcs/Pkt	Pcs/Carton
22.0435	Plastic Type	50	500
22.0452	Disposable type, with 2 fastener batteries	50	500
22.0444	With tongue depressor holder	50	500

Cat.No	Type & Size	Pcs/Pkt	Pcs/Carton
24.0911	Oral	12	720
24.0925	Rectal	12	720
24.0937	Armpite-Medium Size	12	720
24.0939	Armpite-Large Size	12	720

TOURNIQUET

Different types and colours
with safety closure
Single packed



Cat.No	Colour	Pcs/Pkt	Pcs/Carton
22.0524	Blue colour	100	1000
22.0525	Green colour	100	1000
22.0526	Red colour	100	1000
22.0532	Children type	100	1000

SCISSORS & TWEZEERS

Unevirsal Scissors (Bandage Scissors), 19 cm
Scissors with blue color plastic handle,11 cm
Scissors with green color plastic handle, 9 cm
Scissors with red color plastic handle ,13 cm
Plastic Tweezers, Blue Colour.13 cm
Plastic Tweezers ,Sterile ,single PE package, Blue Colour.13 cm



Cat.No	Size	Pcs/Pkt	Pcs/Carton
24.0657	19 cm	24	240
24.0655	11 cm	36	360
24.0651	9 cm	48	480
24.0653	13 cm	36	360
24.0605	13 cm	100	1000
24.9605	13 cm	100	1000

FIRST AID KIT

Pocket Type

Plastic blue colour kit

Small Pocket type recommended for; Cars, Camping, Offices, etc,,,

Dim. Mm 128 x 95 x 45



CONTENTS:

- | | |
|-------------------------------------|------------------------------|
| 1 Adhesive Plaster 1.25 cm x 4.5m | 1 Gauze Bandage 7.5 cm x 5 m |
| 1 Adhesive dressing 6 x 8 cm | 1 Pair of Scissors |
| 5 Wound adhesive plaster 19 x 72 mm | 1 Burn gel 3.5 g |
| 6 Alcohol Swab | |
| 5 Gauze Swab sterile 5 x 5 cm | |

Cat.No	Type	Pkt/Carton
26.0105	Empty	72
26.0106	Filled	72

FIRST AID KIT

Personal Type

First aid Personal Kit

Small type recommended for; Cars, Schools, Offices, etc,,,

Dim.mm 235 x 190 x 43



CONTENTS:

- | | |
|-----------------------------------|------------------------------|
| 1 Adhesive Plaster 1.25 cm x 4.5m | 1 Cotton Balls 5g |
| 1 Adhesive dressing 6x8 cm | 10 Alcohol Swab |
| 2 Gauze bandage 5cmx5m | 1 Sterile Eye Pad |
| 1 Gauze bandage 10cmx5m | 1 Burn gel 3.5 g |
| 5 Gauze Swab sterile 7.5x7.5 cm | 1 Plastic tweezers |
| 10 Wound adhesive plaster 19x72mm | 1 Pair of Scissors |
| 2 Pairs Latex gloves | 1 First Aid Guidance Leaflet |

Cat.No	Colour & Type	Pkt/Carton
26.0127	Orange color-EMPTY.	24
26.0128	Orange color-FILLED.	24
26.0129	Blue color-EMPTY.	24
26.0130	Blue color-FILLED.	24

MEDICAL
DISPOSABLES
PRODUCTS



T H E G E N T L E T O U C H



DISPOSABLE STERILE SYRINGES

Quality disposable syringes, Medical grade raw materials

Latex free piston

Availability:

different sizes, 1,2,3,5,10,20,30,and 50 - 60 ml

Luer Lock /Lour Slip, with /without needle, PE /Blister packing.



	Cat.No	Volume	Gauge	Pcs/Pkt	Pcs/Carton
Luer Lock, Blister packing	31.9201	1 ml	G 27 x 1/2"	100	3000
	31.9219	3 ml	23 G x 1 1/4"	100	2400
	31.9233	5 ml	22 G x 1 1/2"	100	1800
	31.9241	10 ml	21 G x 1 1/2"	100	1200
	31.9251	20 ml	21 G x 1 1/2"	50	400
Luer Lock, PE packing	31.9101	1 ml	G 27 x 1/2"	100	2000
	31.9121	3 ml	23 G x 1 1/4"	100	2000
	31.9135	5 ml	22 G x 1 1/2"	100	1600
	31.9141	10 ml	21 G x 1 1/2"	100	1200
	31.9151	20 ml	21 G x 1 1/2"	100	600

INSULIN SYRINGES

High quality Insulin Syringes for diabetes patients care,

Fine Needle for comfortable use,

Clear graduation with smooth piston movement for accurate dosage,

Singly-Packed, Sterile



Cat.No	Volume	Gauge	Pcs/Pkt	Pcs/Carton
31.9005	1mm/cc	G 29 x 12/	100	2000
	U-100 Insulin			



DISPOSABLE STERILE NEEDLE

Sterile Hypodermic Needle, For single use

100 Pcs/Pkt - 100 Pkt/Case

Sterilized by E.O., Pyrogen Free, Toxic Free, Single blister packing



Cat.No	Gauge	Pcs/Pkt	Pcs/Carton
31.9603	G18 34/100"	100	10,000
31.9606	G18 1 12/100"	100	10,000
31.9616	G19 1 12/100"	100	10,000
31.9624	G20 1"	100	10,000
31.9626	G20 1 12/100"	100	10,000
31.9634	G21 1"	100	10,000
31.9635	G21 1 14/100"	100	10,000
31.9636	G21 1 12/100"	100	10,000
31.9644	G22 1"	100	10,000
31.9645	G22 1 14/100"	100	10,000
31.9646	G22 1 12/100"	100	10,000
31.9654	G23 1"	100	10,000
31.9655	G23 1 14/100"	100	10,000
31.9656	G23 1 12/100"	100	10,000
31.9664	G24 1"	100	10,000
31.9666	G24 1 12/100"	100	10,000
31.9672	G25 58/100"	100	10,000
31.9674	G25 1"	100	10,000
31.9683	G26 34/100"	100	10,000
31.9691	G27 12/100"	100	10,000

I.V CATHETERS

3 Types available:

- I.V Catheter butterfly type
- I.V Catheter with Port & Wings
- I.V Catheter Pen Type

Smooth catheter-cannula connection, secure closure injection port, Single blister packing, Sterile



	Cat.No	Gauge	Colour	Pcs/Pkt	Pcs/Carton
butterfly type	14.9102	14G	orange	50	1000
	14.9103	16G	grey	50	1000
	14.9104	17G	white	50	1000
	14.9105	18G	green	50	1000
	14.9106	20G	pink	50	1000
	14.9107	22G	blue	50	1000
	with Port & Wings	14.9112	14G	orange	50
14.9113		16G	grey	50	1000
14.9114		17G	white	50	1000
14.9115		18G	green	50	1000
14.9116		20G	pink	50	1000
14.9117		22G	blue	50	1000
Pen Type		14.9122	14G	orange	50
	14.9123	16G	grey	50	1000
	14.9124	17G	white	50	1000
	14.9125	18G	green	50	1000
	14.9126	20G	pink	50	1000
	14.9127	22G	blue	50	1000



SCALP VEIN SET

Colour coded,
 Flexible wings for secure placement,
 Stainless steel ultra thin needle,for patient comfort,
 Luer Lock connector with Medical grade PVC tube,
 Singly Packed
 Sterile



Cat.No	Gauge	Colour	Pcs/Pkt	Pcs/Carton
14.9152	18G	pink	50	1000
14.9153	19G	beige	50	1000
14.9154	20G	yellow	50	1000
14.9155	21G	green	50	1000
14.9156	22G	black	50	1000
14.9157	23G	blue	50	1000
14.9158	24G	violet	50	1000
14.9159	25G	orange	50	1000

I.V ADMINISTRATION SET

with wings

20 drops/ml, with or without Injection port, with or without Filter
 Stainless steel ultra thin needle,
 Strong smooth Spike for optimal administration with large roller clamp,
 Singly Packed, Sterile.

With Injection Port (without filter)



With Injection Port (With filter)



Without injection port (without filter)



Cat.No	Type	Pcs/Pkt	Pcs/Carton
14.9202	With Injection Port (without filter)	25	400
14.9204	With Injection Port (With filter)	25	400
14.9212	Without injection port (without filter)	25	400
14.9214	Without injection port (With filter)	25	400

I.V ADMINISTRATION SET

without wings

20 drops/ml
 with or without Injection port
 with or without Filter
 Stainless steel ultra thin needle
 Strong smooth Spike for optimal administration with large roller clamp,
 Singly Packed,
 Sterile.



With Injection Port (With filter)




Without injection port (without filter)

Cat.No	Type	Pcs/Pkt	Pcs/Carton
14.9206	With Injection Port (without filter)	25	400
14.9208	With Injection Port (With filter)	25	400
14.9216	Without injection port (without filter)	25	400
14.9218	Without injection port (With filter)	25	400

BLOOD TRANSFUSION SET

20 drops/ml,
 Integral blood filter,
 Stainless steel ultra thin needle,
 Strong smooth Spike for optimal administration,
 With large roller clamp,
 Singly Packed,
 Sterile.



Cat.No	Type	Pcs/Pkt	Pcs/Carton
14.9411	20 drops/ml	25	400

THREE WAY STOPCOCK

For safe and satable infusion
 Transparent body for clear monitoring
 Arrow indicates the direction of flow,
 360° rotation, with luer lock connectors,
 Singly-Packed, Sterile



Cat.No	Type	Pcs/Pkt	Pcs/Carton
14.9355	Three way stopcock	50	1000
14.9357	Three Way stopcock with extension tube	50	1000
14.9359	Extention Tube	50	1000

SHARP CONTAINERS

Sharps Disposal Containers
 Multi-Purpose Sharp Containers
 Offers a convenient and safe way of disposing contaminated needles,
 syringes, Scalpel Blades,Lancets
 and all types of sharps and medical waste,
 with range of 5 different sizes, and various types of openings.



Cat.No	Type	Pcs/Pkt	Pcs/Carton
32.0302	Round Type-Yellow colour with Rotor opening.	0.5 Ltr.	90
32.0313	Yellow colour-With Rotor opening	2.2 Ltr.	30
32.0315	Yellow colour-With Rotor opening	6.2 Ltr.	24
32.0316	Yellow colour-With Rotor opening,New Type	7 Ltr.	24
32.0317	Yellow colour-Multi purpose with Horizontal drop opening lid	10 Ltr.	12
32.0319	with horizontal drop opening lid.	15 Ltr.	12

DRESSING PACKS

DISPOSABLE DRESSING PACK

- 1 Pec Tray with 3 compartments
- 5 Pcs cotton balls
- 5 Pcs non-woven swabs 30 g, 5x5 cm-4 p
- 4 Pcs -Plastic forceps 11 cm
- 1 Pec Surgical blade 11#
- 1 Pec Alcohol Swab
- 2 Pcs Examination latex gloves, Size M
- 1 Pec Film coated 50 x 50 cm, blue colour
- 1 Pec Sterile Field 60 x 60 cm
- All Packed in Sterilization packing pouch



Cat.No	Size	Pack/Carton
15.9003	Medium	40

DRESSING PACKS

DISPOSABLE DRESSING PACK

- 1 Pec Tray with 3 compartments
- 5 Pcs cotton balls
- 5 Pcs non-woven swabs 30 g, 7.5 x 7.5 cm-4 p
- 4 Pcs -Plastic forceps 11 cm (green & yellow)
- 1 Pec disposable waste bag with adhesive edge
- 1 Pec Film coated 45 x 45 cm, white colour
- 1 Pec Sterile Field 50 x 50 cm
- All Packed in Sterilization P.E bag



Cat.No	Size	Pack/Carton
15.9005	Medium	100

DRESSING PACKS

DISPOSABLE DRESSING PACK

- 1 Pec Tray with 3 compartments
- 8 Pcs cotton balls
- 5 Pcs non-woven swabs 30g, 7.5 x 7.5 cm-4 p
- 3 Pcs -Plastic forceps 11 cm (green & yellow)
- 1 Pec Film coated 45 x 45cm, white colour
- 1 Pec Sterile Field 50 x 50 cm
- All Packed in Sterilization P.E bag



Cat.No	Size	Pack/Carton
15.9007	Medium	100

UMBILICAL CORD CLAMP

- Secure, non-release closure,
- Made of medical grade raw material,
- Sterile,
- Singly-packed in peel-pouch



Cat.No	Colour	Pcs/Pkt	Pcs/Carton
20.9212	White	100	2000

STERILE DISPOSABLE SCALPEL

Surgical scalpels made from high-grade Stainless Steel,
 Plastic handle designed for better feel and performance
 Protective plastic cover,
 Sterile,
 Packed in aluminium foil,
 10 Pcs per dispenser box,



Cat.No	Sizes	Pcs/Pkt
12.9110	10	10
12.9111	11	10
12.9112	12	10
12.9115	15	10
12.9120	20	10
12.9121	21	10
12.9122	22	10
12.9123	23	10
12.9124	24	10
12.9125	25	10

STERILE SURGICAL BLADES

Made from high-grade Stainless Steel,
 Sterile,
 Packed in aluminium foil,
 100 Pcs per dispenser box,



Cat.No	Sizes	Pcs/Pkt
12.9010	10	100
12.9011	11	100
12.9012	12	100
12.9015	15	100
12.9020	20	100
12.9021	21	100
12.9022	22	100
12.9023	23	100
12.9024	24	100
12.9025	25	100

DIPOSABLE UNDERPADS

Absorbent core is made of long virgin pulp fibers for maximum absorbency,
With antibacterial for more protection,
Non-woven top sheet for optimum dryness and softness,
Polyethylene waterproof back sheet for leakage protection,



Cat.No	Sizes	Pcs/Pkt	Pcs/Carton
35.0108	60x40 cm	30	180
35.0110	60x60 cm	30	180
35.0112	60x90 cm	30	180

PAEDIATRIC URINE BAG COLLECTOR

Capacity: 100 ml,
with adhesive area,
for short term urine collection,
Suitable for girls and boys,
Sterile,
Singly-packed in peel-pouch



Cat.No	Size	Pcs/Pkt
19.9052	100 ml	100

SHOE COVER

CPE disposable Shoe Covers,
Ideal for Operating theatre, Laboratory
and food industry
Blue Color,
Packed in P.E bag
100 PCS/Pkt



Cat.No	Type	Pcs/Pkt	Pcs/Carton
18.0503	2.3 g, 36 X 15 cm	100	2000

I.D BRACELET

Identity Band
PVC material
Available colours;
White,Blue, Pink, Red and Yellow



	Cat.No	Colour	Pcs/Pkt	Pcs/Carton
Adult	18.0612	Blue	100	5000
	18.0613	Light Blue	100	5000
	18.0614	white	100	5000
	18.0615	Red	100	5000
	18.0616	Yellow	100	5000
	18.0617	Pink	100	5000
	Children	18.0622	Blue	100
18.0623		Light Blue	100	5000
18.0624		white	100	5000
18.0625		Red	100	5000
18.0626		Yellow	100	5000
18.0627		Pink	100	5000

CAST PADDING

Orthopaedic Padding
 Polyester Type
 90 g /sq.m
 White colour



Cat.No	Size	Pcs/Pkt	Pcs/Carton
18.0710	5 cm x 270 cm	12	600
18.0716	7.5 cm x 270 cm	12	480
18.0722	10 cm x 270 cm	12	360
18.0728	15 cm x 270 cm	12	240
18.0734	20 cm x 270 m	12	180

PLASTER OF PARIS

P.O.P bandage
 450 g /sq.m
 Drying time; 2-8 minutes



Cat.No	Size	Pcs/Pkt	Pcs/Carton
18.0802	2"x 270 cm	12	240
18.0803	3"x 270 cm	12	240
18.0804	4"x 270 cm	12	120
18.0806	6"x 270 cm	6	120
18.0808	8"x 270 cm	6	60
18.0810	10"x 270 cm	6	48
18.0812	12"x 270 cm	6	36

SURGICAL MASK

Ear Loop

Surgical Disposable Face Masks

High breathability

Nose adapter for correct and comfort face adjustment

High Filter action (BFE ≥ 95 %)

Ear Loop Type, 50 Pcs/Pkt



Cat.No	Type	Pcs/Pkt	Pcs/Carton
16.0206	2 Ply, ear loop type, blue colour	50	2000
16.0207	2 Ply, ear loop type, green colour	50	2000
16.0208	2 Ply, ear loop type, White colour	50	2000
16.0210	3 Ply, ear loop type, blue colour	50	2000
16.0211	3 Ply, ear loop type, green colour	50	2000
16.0212	3 Ply, ear loop type, White colour	50	2000
16.0280	3 Ply, ear loop, blue colour Excellent Filtration (BFE 99 %)	50	2000

SURGICAL MASK

Tie Type


Surgical Disposable Face Masks

High breathability

Nose adapter for correct and comfort face adjustment

High Filter action (BFE ≥ 95 %)

Tie Type, 50 Pcs/Pkt



Cat.No	Type	Pcs/Pkt	Pcs/Carton
16.0252	2 Ply tie type, blue colour	50	2000
16.0253	2 Ply tie type, green colour	50	2000
16.0254	2 Ply tie type, white colour	50	2000
16.0260	3 Ply tie type, blue colour	50	2000
16.0261	3 Ply tie type, green colour	50	2000
16.0262	3 Ply tie type, white colour	50	2000
16.0282	3 Ply tie type, blue colour Excellent Filtration (BFE 99 %)	50	2000

SURGEON CAPS

Surgical Disposable Head Cap
 SMS Material
 100 Pcs/ Dispenser Pack



	Cat.No	Type	Pcs/Pkt	Pcs/Carton
Blue colour	16.0305	SMS 20 gm - Tie Type	100	1000
	16.0307	SMS 20 gm - Elastic Type	100	1000
	16.0309	SMS 40 gm - Tie Type	100	1000
	16.0311	SMS 40 gm - Elastic Type	100	1000
Green Colour	16.0315	SMS 20 gm - Tie Type	100	1000
	16.0317	SMS 20 gm - Elastic Type	100	1000
	16.0319	SMS 40 gm - Tie Type	100	1000
	16.0321	SMS 40 gm - Elastic Type	100	1000

SURGEON CAPS

Surgical Disposable Head Cap
 Non-woven Material
 Light and breathable
 100 Pcs/ Dispenser Pack



	Cat.No	Type	Pcs/Pkt	Pcs/Carton
Blue colour	16.0352	Non woven, Tie Type, 30 gm.	100	1000
	16.0354	Non woven, Elastic Type, 30 gm.	100	1000
Green Colour	16.0366	Non woven, Tie Type, 30 gm.	100	1000
	16.0368	Non woven, Elastic Type, 30 gm.	100	1000



NURSE HEAD CAP

Bouffant Type

Bouffant Nurse Cap
 Non-woven material
 Light and breathable
 With rubber head-band
 100 Pcs/ Dispenser Pack



Cat.No	Type	Pcs/Pkt	Pcs/Carton
16.0463	20 g , Blue colour , 21"	100	1000
16.0465	20 g , Green colour , 21"	100	1000
16.0467	20 g , White colour , 21"	100	1000

CLIP CAP

Mushroom Type

Head Cap, for nurses, visitors and for industrial use,
 Non-woven material
 Light and breathable
 Universal size with double elastic bands for extra protection
 100 Pcs/Pkt



Cat.No	Type	Pcs/Pkt	Pcs/Carton
16.0411	14 g, White Colour , 21"	100	2000
16.0413	14 g, Blue Colour , 21"	100	2000
16.0415	14 g, Green Colour , 21"	100	2000

STERILE GAUZE SWAB

Made of 100 % cotton

Soft and Durable with high level of absorbency and low linting,

Different meshes, sizes, plies, and packing are available



	Cat.No	Size	Pcs/Pkt	Pcs/Carton
40s/19x15 1 Pcs/Blister	06.9201	5 x 5 cm - 8 Plies	100	2000
	06.9202	7.5 x 7.5 cm - 8 Plies	100	2000
	06.9203	10 x 10 cm - 8 Plies	100	2000
	06.9401	5 x 5 cm - 12 Plies	100	1600
	06.9402	7.5 x 7.5 cm - 12 Plies	100	1600
	06.9403	10 x 10 cm - 12 Plies	100	1600
40s/19x15 5 Pcs/Blister	06.9231	5 x 5 cm - 8 Plies	100	2000
	06.9232	7.5 x 7.5 cm - 8 Plies	100	2000
	06.9233	10 x 10 cm - 8 Plies	100	2000
	06.9431	5 x 5 cm - 12 Plies	100	1600
	06.9432	7.5 x 7.5 cm - 12 Plies	100	1600
	06.9433	10 x 10 cm - 12 Plies	100	1600

STERILE GAUZE SWAB

with X-Ray Detectable

Made of 100 % cotton

with X-Ray thread

Soft and Durable with high level of absorbency and low linting,

Different meshes, sizes, plies, and packing are available



	Cat.No	Size	Pcs/Pkt	Pcs/Carton
40s/19x15 1 Pcs/Blister	07.9201	5 x 5 cm - 8 Plies	100	2000
	07.9202	7.5 x 7.5 cm - 8 Plies	100	2000
	07.9203	10 x 10 cm - 8 Plies	100	2000
	07.9401	5 x 5 cm - 12 Plies	100	1600
	07.9402	7.5 x 7.5 cm - 12 Plies	100	1600
	07.9403	10 x 10 cm - 12 Plies	100	1600
40s/19x15 5 Pcs/Blister	07.9231	5 x 5 cm - 8 Plies	100	2000
	07.9232	7.5 x 7.5 cm - 8 Plies	100	2000
	07.9233	10 x 10 cm - 8 Plies	100	2000
	07.9431	5 x 5 cm - 12 Plies	100	1600
	07.9432	7.5 x 7.5 cm - 12 Plies	100	1600
	07.9433	10 x 10 cm - 12 Plies	100	1600

LAPAROTOMY SPONGES

Lap Sponges

Made of 100 % cotton

High Level of absorbency and durability, Chlorine free,

5 Pcs/Blister Packing

available as:

Sterile and non sterile, Pre-Washed and non-washed

with or without Loop, with or without X-Ray thread

Different meshes, Sizes & Plies, Sterile, with Loop and X-Ray Detectable



	Cat.No	Size	Pcs/blister	Pcs/Carton
8 Plies	08.9409	22.5 x 22.5 cm	5	250
	08.9411	30 x 30 cm	5	250
	08.9413	40 x 40 cm	5	250
	08.9415	45 x 45 cm	5	250
12 Plies	08.9609	22.5 x 22.5 cm	5	250
	08.9611	30 x 30 cm	5	250
	08.9613	40 x 40 cm	5	250
	08.9615	45 x 45 cm	5	250

SURGICAL GOWNS

Sterile Surgical Gowns

Blue Colour, Double Packing.

Made of SMS material for maximum protection

water, Oil, Alcohol, and blood repellent,

Comfortable, Soft and light weight for long procedures,

Inner Packing, 1 Pec/Pouch

Outer Packing, 50 Pcs/Ctn



Cat.No	Dim	Size	Pcs/Carton
04.9204	115 x 130 cm	Medium	50
04.9205	120 x 150 cm	Large	50
04.9206	135 x 160 cm	X-Large	50
04.9207	135 x 200 cm	XX-Large	50
04.9208	135 x 230 cm	XXX-Large	50

ISOLATION GOWNS

Non Sterile Isolation Gown
 Non-Woven Material
 Elastic cuff
 Uni-Size
 Inner Packing: 10 Pcs/Bag
 Outer Packing: 10 bags/Ctn



HMS Surgical Instruments

Cat.No	Colour	Size	Pcs/Pkt	Pcs/Carton
04.0653	Blue	UNI-Size	10	100
04.0663	Yellow	UNI-Size	10	100



NOTES

A series of horizontal dotted lines for taking notes, spanning the width of the page.



Halley Medical Supplies Ltd.

1406 Coventry Road, South Yardley,
 Birmingham B25 8AE UK.
 Tel: + 44 121 7650214
 Fax: + 44 121 7650215
 sales@halleymedical-uk.ae
 www.halleymedical.co.uk



NOTES

A series of horizontal dashed lines for writing notes.